

SFWRENG 3K04

Deliverable 1

Group 28

Fall 2025

Contents

1 Group Members	2
2 Part 1	2
2.1 Introduction	2
2.2 Requirements	3
2.3 Design	4
3 Part 2	15
3.1 Requirements Potential Changes	15
3.2 Design Decision Potential Changes	16
3.3 Module Description	17
3.4 Testing	19
3.5 GenAI Usage	21
4 General Notes	4

1 Group Members

Aaditya Anil, anila4, 400520758

Deev Patel, pated201

Shaun Plassery, plasers,

Dylan Manamendra, manamenk, 400509187

2 Part 1

2.1 Introduction

Arrhythmia is a conduction disorder in which the heart's electrical signal is disrupted inducing, irregular heart rhythm or rate. Individuals with a form of arrhythmia are at

great risk of heart attacks, cardiac arrest, organ failure, and stroke (NHLBI, NIH – March 2022). To counteract this, a pacemaker is used to detect when the heart rate of the user drops, then administer a controlled electrical signal to the muscles surrounding the heart to restabilize its natural rhythm. The pacemaker being designed currently is designed to treat a form of arrhythmia known as Bradycardia, in which the heart rhythm drops to a dangerously low rate (<60 bpm). The pacemaker is equipped with a suite of sensors to detect information on the user's heart rate, checking if said info is in line with the conditions of those experiencing an arrhythmic episode. Once an arrhythmic episode is detected, the pacemaker administers a controlled electrical signal to the surrounding muscular tissue, forcing the heart to beat at an appropriate rate.

The Device Control Module (DCM) is the software interface that enables clinicians to configure and monitor the pacemaker device. It provides user authentication, pacing-mode configuration, programmable parameter entry, and placeholder functionality for egram visualization. For this deliverable, only the presentation layer (front-end) of the DCM has been implemented; communications between the DCM and pacemaker hardware are planned to be incorporated in the next phase

This first deliverable is to document the initial requirements of the project, and ensure that the group members not only understand the function of the pacemaker, but are aware of the necessary specifications needed to design and construct the state flow diagrams of the pacemaker and the graphical user interface (GUI) for the DCM.

2.2 Requirements

- Pulse width and refractory period timing of atrium and ventricle
 - Information on the low-rate limit and upper-rate limit (LRL, URL)
 - Hysteresis (irregularity in pulse rate) state detection
 - Pacing capacitor with sufficient power source
- o AOO: grounding control (Pace Cap/Ventricle/Atrium), pace control (Atrium, Ventricle), DCM control, LRL, URL, atrial pulse width
 - o VOO: grounding control, pace control (Atrium, Ventricle), DCM control, LRL, URL, ventricle pulse width
 - o AAI: grounding control, pace control (Atrium, Ventricle), DCM control, LRL, URL, hysteresis state, atrium pulse width, atrial refractory period

- o VVIL grounding control, pace control (Atrium, Ventricle), DCM control, LRL, ventricle pulse width, hysteresis state, ventricle refractory period

The following summarizes how the implemented DCM fulfills **section 3.2.2 of the PACEMAKER System Specification**:

Req ID	Description	Type	Trace to Design Section
R1	The DCM shall allow registration and login for up to 10 users.	Functional	DCM Design § User Utilities
R2	The DCM shall display and edit programmable parameters (LRL, URL, AA, APW, VA, VPW, VRP, ARP).	Functional	DCM Design § Parameter Manager
R3	The DCM shall provide a status indicator for device connection.	Functional	DCM Design § DeviceManager
R4	The DCM shall save all parameter data locally in JSON format.	Functional	DCM Design § File Persistence
R5	The DCM shall indicate when a different device is detected.	Functional	DCM Design § DeviceManager
R6	The system shall validate all numeric parameters against range limits.	Non-Functional (Safety)	DCM Design § Parameter Validation

Additional derived requirements from the assignment:

- Store ≤ 10 local users.
- Provide input validation for parameter ranges (Section PARAM_RANGES).
- Allow saving and loading of parameter files.
- Maintain a placeholder for future Egram recording and plotting.

2.3 Design

The state flow model, designed using Simulink, illustrates the direction and behavior of data within the pacemaker device, highlighting which states use specific information (LRL, atrium pulse width, etc.) for their operation.

The pacemaker is fed information regarding the LRL and URL, the mode of operation, as well as sensor data detecting the atrial and ventricle compression and pulse width, atrial and ventricle refractory periods, and whether the heart rate has entered hysteresis.




Figure 1: Pacemaker State Machine with external input data

The mode variable is responsible for setting the state of the pacemaker, and subsequently the means in which it functions.




Figure 2: Mode value used as state control for pacemaker

Within each state is the parameters and data used for the specific function, as well as transition conditions between specific operations within the state. When entering the AOO state, the state flow starts a timer; upon reaching the computed LRL interval it declares atrial pace control, ATR_PACE_CTRL, for Atr_PulseWidth, then resets the interval. Additionally, it allows for the pacing capacitor to be charged up. No sensing is performed (inhibit = none), and the pacemaker enters “no response” behavior. AOO applies only to the atrial activity.




Figure 3: AOO state of pacemaker

When entering the VOO state, the state flow once again starts a timer; upon reaching the computed LRL interval it declares the ventricle pace control, VENT_PACE_CTRL, for Vent_PulseWidth, then resets the interval whilst also charging the pacing capacitor. The function of VOO is similar to AOO, except applying to the ventricle activity as opposed to the atrial activity.




Figure 4: VOO state of pacemaker

When entering the AAI state, the controller starts the LRL timer. If ATR_CMP_DETECT (Boolean to detect atrial compression) occurs before the timer elapses and not in refractory, the timer is reset, and no pace is delivered. If the timer completes without an atrial sense, the state flow declares ATR_PACE_CTRL for Atr_PulseWidth. The pacemaker will also discharge the pacing capacitor, and allow it to charge up again. This once again, is only with regards to atrial activity




Figure 5: AAI state for pacemaker

When entering the VVI state, the pacemaker computes $f_{waitInterval}$ based on $p_{lowrateInterval}$ (found using LRL) plus hysteresis if enabled, account for ventricular refractory periods (VRP), and detect spontaneous ventricular sense in the wait interval. If detected, a ventricular pace is delivered (VENT_PACE_CTRL) by pacing capacitors, so long as neither sensed nor paced events occur within the wait interval. The chart implements In_vPace , In_pVRP , and the c_vp output logic as described in the pacemaker document.




Figure 6: VVI state of the pacemaker

In all states of the pacemaker, conditions such as `after()` and values such as LRL, URL, vent/atr pulsewidth, and internal timers are used to implement intervals and pulse widths (e.g., `after(LRL - Vent_Pulse_Width, msec)` for pacing decision timing and `after(p_vPaceWidth, msec)` to clear a pacing pulse). These match the timings and tolerances for pulse width and refractory intervals labeled on the pacemaker documentation. Additionally, expectations such as asynchronous operation of AOO and VOO states, AAI and VVI control, manual control via the DCM, and Hardware hiding to correct board ports (for simulation of an actual heart) have been met in accordance to the deliverable 1 expectations.




Figure 7:Hardware hiding of outputs to correct ports in microcontroller




Figure 8: Port connection of specific waveforms within the microcontroller

The DCM follows a **modular, object-oriented structure**, divided into self-contained components for user management, device management, parameter storage, and data visualization.

Major Modules

- **WelcomeApp (class)** – Main Tkinter application controller; manages frame navigation and state.
- **DeviceManager (class)** – Tracks the most recently interrogated device to support requirement (7).
- **User Management Utilities** – Handle registration and authentication with password hashing and salt.
- **Parameter Management Subsystem** – Validates and stores programmable parameters per mode.
- **Egram Subsystem** – Defines JSON schema and save/list/preview functions for recorded Egrams.

Below the User Login → Main DCM → Parameter Manager → Egram Module block diagram can be found.




Figure 9: Block Diagram of DCM

The welcome screen offers three options — Login, Register, or Exit. Registration supports up to 10 users stored locally in `users.json`. Passwords are hashed with SHA-256 and per-user salts for security.

```

385 class WelcomeApp(tk.Tk):
386     # -----
387     # Welcome / login / register
388     # -----
389     def build_welcome(self):
390         f = self.frame_welcome
391         tk.Label(f, text="Welcome!", font=("Arial", 18)).pack(pady=(0, 10))
392         tk.Button(f, text="Login", width=20, command=lambda: self.show_frame(self.frame_login)).pack(pady=5)
393         tk.Button(f, text="Register", width=20, command=lambda: self.show_frame(self.frame_register)).pack(pady=5)
394         tk.Button(f, text="Exit", width=20, command=self.quit).pack(pady=(10, 0))
395
396     def build_login(self):
397         f = self.frame_login
398         tk.Label(f, text="Login", font=("Arial", 14)).grid(row=0, column=0, colspan=2, pady=(0, 8))
399         tk.Label(f, text="Username:").grid(row=1, column=0, sticky="e")
400         self.login_user = tk.Entry(f); self.login_user.grid(row=1, column=1)
401         tk.Label(f, text="Password:").grid(row=2, column=0, sticky="e")
402         self.login_pass = tk.Entry(f, show="*"); self.login_pass.grid(row=2, column=1)
403         self.login_msg = tk.Label(f, text="", fg="red"); self.login_msg.grid(row=3, column=0, colspan=2)
404         tk.Button(f, text="Submit", width=12, command=self.attempt_login).grid(row=4, column=0, pady=8)
405         tk.Button(f, text="Back", width=12, command=lambda: self.show_frame(self.frame_welcome)).grid(row=4, column=1)
406
407     def build_register(self):
408         f = self.frame_register
409         tk.Label(f, text="Register", font=("Arial", 14)).grid(row=0, column=0, colspan=2, pady=(0, 8))
410         tk.Label(f, text="Username:").grid(row=1, column=0, sticky="e")
411         self.reg_user = tk.Entry(f); self.reg_user.grid(row=1, column=1)
412         tk.Label(f, text="Password:").grid(row=2, column=0, sticky="e")
413         self.reg_pass = tk.Entry(f, show="*"); self.reg_pass.grid(row=2, column=1)
414         tk.Label(f, text="Confirm:").grid(row=3, column=0, sticky="e")
415         self.reg_pass2 = tk.Entry(f, show="*"); self.reg_pass2.grid(row=3, column=1)
416         self.reg_msg = tk.Label(f, text="", fg="red"); self.reg_msg.grid(row=4, column=0, colspan=2)
417         tk.Button(f, text="Create", width=12, command=self.attempt_register).grid(row=5, column=0, pady=8)
418         tk.Button(f, text="Back", width=12, command=lambda: self.show_frame(self.frame_welcome)).grid(row=5, column=1)
419
420     def attempt_register(self):

```

Figure 10: Code for login window




Figure 11: Login window

Upon login, the user is brought to the main DCM interface, which contains three panes:

- **Left Pane:** Connection status indicator, device dropdown, and refresh/disconnect buttons.

- **Middle Pane:** Dynamic parameter fields for all programmable parameters.
- **Right Pane:** Communication log panel displaying DCM events.

The top bar shows the logged-in user and a logout button.

```

389 class WelcomeApp(tk.Tk):
499     def _build_dcm(self):
541         # Right: packet log
542         right = tk.Frame(f)
543         right.pack(side="left", padx=(12,0))
544         tk.Label(right, text="Comm Log", font=("Arial", 10)).pack(pady=(0,4))
545         self.log_text = tk.Text(right, width=48, height=16, state="disabled")
546         self.log_text.pack(pady=(4,0))
547
548         # placeholder for egram viewer + quick egram actions
549         frame_eg = tk.Frame(f)
550         frame_eg.pack(side="bottom", pady=(12,0))
551         tk.Label(frame_eg, text="Egram Viewer (placeholder)", fg="gray").pack(side="left")
552         # quick egram snapshot / list buttons
553         try:
554             tk.Button(frame_eg, text="Save Egram Snapshot", command=self._ui_save_gram).pack(side="left", padx=(8,4))
555         except Exception:
556             # if ttk isn't imported, fall back
557             btn = tk.Button(frame_eg, text="Save Egram Snapshot", command=self._ui_save_gram)
558             btn.pack(side="left", padx=(8,4))
559         try:
560             tk.Button(frame_eg, text="List Egrams", command=self._ui_list_egms).pack(side="left", padx=(4,0))
561         except Exception:
562             btn2 = tk.Button(frame_eg, text="List Egrams", command=self._ui_list_egms)
563             btn2.pack(side="left", padx=(4,0))
564
565         # initial device list
566         self._refresh_device_list()
567         # trace selection
568         self.device_var.trace_add("write", self._on_device_selected)
569
570

```

Figure 12: Code for Main DCM Interface




Figure 13: Main Page for DCM

Each mode (e.g., AOO, VVI, DDD) enables only its relevant parameters based on the mapping in PARAMS_BY_MODE.

Entries are validated through `_clamp_to_range()` and `_is_number()` before saving.

If a parameter is out of range (e.g., LRL > URL or amplitude > 7.5 V), an error dialog appears via `messagebox.showerror()`.

```

385 class WelcomeApp(tk.Tk):
499     def build_dcm(self):
526         # Middle: parameter viewer / editor
527         mid = tk.Frame(f)
528         mid.pack(side="left")
529         tk.Label(mid, text="Programmable Parameters", font=("Arial", 12)).pack()
530
531         self.params_frame = tk.Frame(mid, bd=1, relief="sunken", padx=6, pady=6)
532         self.params_frame.pack(pady=(6,4))
533
534         # Parameter entries will be created dynamically
535         self.param_vars = {}
536         self.param_entries = {}
537
538         btns = tk.Frame(mid); btns.pack(pady=6)
539         tk.Button(btns, text="Save to File", command=self._save_params_to_file).pack(side="left", padx=6)
540         tk.Button(btns, text="Show Last Device Info", command=self._show_last_device_info).pack(side="left")
541
542         # Right: packet log
543         right = tk.Frame(f)
544         right.pack(side="left", padx=(12,0))
545         tk.Label(right, text="Comm Log", font=("Arial", 10)).pack(pady=(0,4))
546         self.log_text = tk.Text(right, width=48, height=16, state="disabled")
547         self.log_text.pack(pady=(4,0))
548
549         # placeholder for egram viewer + quick egram actions
550         frame_eg = tk.Frame(f)
551         frame_eg.pack(side="bottom", pady=(12,0))
552         tk.Label(frame_eg, text="Egram Viewer (placeholder)", fg="gray").pack(side="left")
553         # quick egram snapshot / list buttons

```

Figure 14: Code for Parameter Editor




Figure 15: Parameter Editor

A status label in the left pane changes color to represent connection state:

- Gray = Disconnected
 - Green = Connected
- If a new device is selected, the program compares the device ID and prompts the user to confirm switching.

```
385 class WelcomeApp(tk.Tk):
499     def build_dcm(self):
511         tk.Label(left, text="Connection").pack()
512         self.status_label = tk.Label(left, text="Disconnected", bg="lightgray", width=18)
513         self.status_label.pack(pady=(4,8))
514
515         # Device selector dropdown
516         tk.Label(left, text="Choose device file:").pack(pady=(6,2))
517         self.device_var = tk.StringVar(value="(none)")
518         self.device_menu = tk.OptionMenu(left, self.device_var, "(none)")
519         self.device_menu.config(width=22)
520         self.device_menu.pack()
521         tk.Button(left, text="Refresh Device List", width=18, command=self._refresh_device_list).pack(pady=(6,2))
522         # Disconnect placed under the device selector as requested
523         tk.Button(left, text="Disconnect", width=18, command=self._dcm_disconnect).pack(pady=(6,2))
524
525         # Middle: parameter viewer / editor
526         mid = tk.Frame(f)
527         mid.pack(side="left")
528
```

Figure 16: Code for status label




Figure 17: Status Label

The current version implements the Egram data structure and simulated waveform generation. Recorded data are saved as JSON files and previewed using Tkinter Canvas. Future iterations will support live signal streaming and real-time plotting.

```

385 class WelcomeApp(tk.Tk):
499     def _build_dcm(self):
548         # placeholder for egram viewer + quick egram actions
549         frame_eg = tk.Frame(f)
550         frame_eg.pack(side="bottom", pady=(12,0))
551         tk.Label(frame_eg, text="Egram Viewer (placeholder)", fg="gray").pack(side="left")
552         # quick egram snapshot / list buttons
553         try:
554             tk.Button(frame_eg, text="Save Egram Snapshot", command=self._ui_save_egram).pack(side="left", padx=(8,4))
555         except Exception:
556             # if ttk isn't imported, fall back
557             btn = tk.Button(frame_eg, text="Save Egram Snapshot", command=self._ui_save_egram)
558             btn.pack(side="left", padx=(8,4))
559         try:
560             tk.Button(frame_eg, text="List Egrams", command=self._ui_list_egms).pack(side="left", padx=(4,0))
561         except Exception:
562             btn2 = tk.Button(frame_eg, text="List Egrams", command=self._ui_list_egms)
563             btn2.pack(side="left", padx=(4,0))
564

```

Figure 18: Current Code for Egram Viewer




Figure 19: Current Egram Viewer

3 Part 2

3.1 Requirements Potential Changes

Scope reminder (D1): Presentation layer only (Tkinter GUI), with local user auth, local device param files, parameter validation by mode, “different device” detection, and offline egram placeholder.

Likely changes for Deliverable 2 (per assignment brief):

- **R-01 Modes & Parameters (expand):** Extend from minimal modes/params to **all required modes** and the **complete parameter set** (PACEMAKER Table 7; srsVVI §3.1 minimal set).
- **R-02 Communication (new):** Add **serial communication** to transmit/receive parameters and egram data (configure port, baud, timeouts).
- **R-03 Parameter Verification (new):** After sending parameters, **read back** from the pacemaker and verify exact match (include data typing and unit conversions; show mismatch).
- **R-04 Parameter Constraints (update):** Adopt **new ranges and increments**:
 - A/V Amplitude: Off, **0.1–5.0 V @ 0.1 V** step (Nominal 5 V $\pm 12\%$)
 - A/V Pulse Width: **1–30 ms @ 1 ms** step
 - A/V Sensitivity: **0–5 V @ 0.1 V** step ($\pm 2\%$)
- **R-05 Egram Streaming (upgrade):** Replace simulated snapshots with **live egram display** received over serial (A, V, or both channels). Allow pause/resume and snapshots.
- **R-06 Telemetry states (refine):** Extend current “Connected/Disconnected” to **timeout/error/noise/out-of-range** conditions if device exposes these.
- **R-07 Persistence (extend):** In addition to local JSON, maintain **TX/RX verification logs** (timestamped records of sent vs. read-back values).
- **R-08 Safety/UX (add):** Non-blocking UI during serial I/O, clear operator prompts on parameter step violations and verification failures.

3.2 Design Decision Potential Changes

- **DD-01** **Parameter** **Metadata** **Source:**
Current: Optional param_metadata.json overrides in _load_param_metadata().
Change: **Adopt metadata as authoritative** for all ranges, units, and **step sizes**, keeping code defaults as fallback only. Add field for step and display units in labels.
- **DD-02** **Validation** **Model:**
Current: Range checks via _clamp_to_range, numeric checks with _is_number, and **mode-gated** enable/disable via _mode_allowed_params.
Change: Add **step-size enforcement** (e.g., multiples of 0.1 V / 1 ms). Validate **LRL \leq URL** (already present) and any **inter-param constraints** introduced by Table 7.
- **DD-03** **Architecture** **for** **I/O:**
Current: Single-threaded Tkinter app; file-based params and simulated

egram.

Change: Introduce a **SerialLink module** that runs on a **worker thread** (or uses Tk `after()` with non-blocking reads). Use a **thread-safe queue** to post events (RX frames, errors) to the UI.

- **DD-04 Protocol Framing & Types:**
Current: N/A.
Change: Define **wire protocol** (message headers, payloads, checksums/CRC, endianness, signedness, scaling factors). Decide **float vs fixed-point** for amplitude/sensitivity; document mapping.
- **DD-05 Plotting:**
Current: Tk Canvas for static preview from JSON.
Change: Keep Tk Canvas for simplicity **or** switch to `matplotlib` (embed in Tk) for smooth **live streaming**. Ensure **no UI freeze** and bounded memory.
- **DD-06 Storage:**
Current: `users.json`, `*_params.json`, `egrams/*.json`, `last_device.json`.
Change: Add `verifications/*.json` (sent vs readback records), **optional SQLite** later if needed.
- **DD-07 Error Handling & UX:**
Current: Message boxes + status label.
Change: Centralize errors in a **Comm Log** with levels (INFO/WARN/ERROR). Add **retry** and **Reconnect** controls; keep visual badges in status bar.

3.3 Module Description

Pacemaker State Machine SubSystem

- **Purpose:** Implements permanent-state pacing behavior across AOO, VOO, AAI, and VVI states.
- **Key Functions:**
 - Computes pacing intervals based on LRL and URL.
 - Executes asynchronous pacing for AOO and VOO.
 - Implements inhibited pacing using atrial or ventricular sensing in AAI and VVI.
 - Generates pacing pulses of programmable width and duration.

- **Internal State Variables:** Timing counters, pulse-width timers, flags indicating sense detection.
- **Interactions:** Consumes programmable parameters and sensing inputs; outputs pacing control signals.

Microcontroller Parts SubSystem (Hardware Hiding Layer)

- **Purpose:** Translates abstract pacing and sensing control signals into specific GPIO pin assignments.
- **Key Functions:**
 - Maps ATR_PACE_CTRL, VENT_PACE_CTRL to hardware pin blocks.
 - Routes ATR_CMP_DETECT, VENT_CMP_DETECT from hardware front-end to the controller.
- **State/Global Variables:** None; behaves as a combinational interface.
- **Interactions:** Directly interfaces with board digital I/O; isolates Stateflow from hardware changes.

DCM

- **Purpose:** Provides GUI for doctors in OR and user outside medical setting to manually control state of pacemaker and relay reports of user cardiac data to medical professionals when needed
- **Key Functions:**
 - Allows users to login to their account and view patient specific data such as medical and cardiac reports, active status of device and , and dates of previous medical appointments or emergencies.
 - Select specific modes of operation/states the pacemaker should be in at any given time (ex. AOO, VOO, AAI, VVI).
 - Can modify constants of operation within the state flow design, such as LRL and URL, expected pulse width and amplitudes
- **State/Global Variables:** None currently; will manage stored parameter presets or user session configurations in later deliverables.

- **Interactions:** Will interface with the Pacemaker State Machine to supply updated parameters. May also receive telemetry output from the State flow chart for display or logging.

3.4 Testing

Document test cases for each module. Each test case should include:

1. Purpose of the test
2. Input conditions
3. Expected output
4. Actual output
5. Result (Pass/Fail)

Simulink State Flow Testing

Purpose	Input	Expected Output
Test AOO pacing	Mode AOO, no inputs	ATR_PACE_CTRL outputs pulses at correct timing, board light flashes accordingly
Test VOO pacing	Model VOO, no inputs	VENT_PACE_CTRL outputs pulses at correct timing, board light flashes accordingly
Test AAI	Model AAI, atrial sensing (pulse, period)	No atrial pace is triggered
Test VVI	Model VVI, ventricle sensing (pulse, period)	No ventricular pace is triggered

DCM Testing

ID	Purpose	Input / Setup	Expected Output	Actual Output	Result
T-01	User registration and authentication	Register a new user "test1"; attempt duplicate registration and login with correct/incorrect passwords	New user created successfully; duplicate rejected; correct password logs in; wrong password shows red error	Works as expected — login transitions to DCM, duplicates blocked	Pass
T-02	Parameter range validation and save	Load sample_parameters.json, edit LRL = 20 (below 30 bpm) and URL = 300 (above 175 bpm)	Error dialogs appear: "Out of range [30–170 bpm]"; parameters not saved	Both validation messages displayed; no file overwritten	Pass
T-03	Device detection and status feedback	Load device_A_parameters.json, then select device_B_parameters.json	Pop-up warns: "Different pacemaker detected"; status changes to Connected (green) / Disconnected (gray) appropriately	Prompt and color indicators function correctly	Pass
T-04	Egram snapshot and preview	Click "Save Egram Snapshot", then "List Egrams → Preview Selected"	New file appears in /egram; preview shows two traces (blue = Atrial, red = Ventricular)	Snapshot JSON created; preview renders waveform correctly	Pass
T-05	Parameter file save and	Edit URL = 130 bpm; click "Save to File"; reopen DCM and reload file	Updated JSON reflects new value; Comm Log: "Saved parameters to	File saved with correct data; log	Pass

persistence

device_A_params.json
"

and
popup
confirmed

3.5 GenAI Usage

GenAI was used for the formatting of this document